

*Lokum “L-Quest Mini“ Guide Wire*



- “L-Quest“ Extra Stiff Guide Wire **MINI** Where High Stability Is Needed
- High PTFE Coating Performance
- Extra Steerable And Pushable For Increased Control

**New L-Quest Mini in 0.025“ !!!**

# AndraTec Advanced Interventional Technologies



*Lokum*  
L-Quest Guide Wire

- The LOKUM L-Quest MINI Guide Wire is especially designed for extra strong support and high pushability in difficult procedures.
- The high performance PTFE coating combined with extra stiff stainless steel LOKUM L-Quest MINI is an excellent choice to treat difficult vascular and non-vascular lesions.
- The Lokum L-Quest guide wire MINI is intended to use for complex interventional procedures where increased stability and low surface friction of the guidewire is needed.

## Technical Information:

Guide Wire Material:	Extra Tempered Stainless Steel
Coating:	High Performance PTFE
Diameter:	0.025" / 0.635 mm
Distal Tip:	Straight / 3mm J-Tip / 6mm J-Tip
Tip Configuration:	5 cm / 10cm / 15cm Flexible Tip
Shaft Length:	150cm / 260 cm

\* a choice of various sizes will assist your patient's needs.

LOKUM "L-Quest MINI" Order Codes (5 per Box)			
0.025" / 0.635mm	Straight Tip	3mm J-Tip	6mm J-Tip
Tip / Shaft Length			
5 cm / 150 cm	LUMI-05-150-ST	LUMI-05-150-3T	LUMI-05-150-6T
10 cm / 150 cm	LUMI-10-150-ST	LUMI-10-150-3T	LUMI-10-150-6T
15 cm / 150 cm	LUMI-15-150-ST	LUMI-15-150-3T	LUMI-15-150-6T
5 cm / 260 cm	LUMI-05-260-ST	LUMI-05-260-3T	LUMI-05-260-6T
10 cm / 260 cm	LUMI-10-260-ST	LUMI-10-260-3T	LUMI-10-260-6T
15 cm / 260 cm	LUMI-15-260-ST	LUMI-15-260-3T	LUMI-15-260-6T



Distributed By:



AndraTec GmbH  
Simmerner Str. 70  
D-56075 Koblenz  
Tel: +49 261 29675 172  
Fax: +49 261 29675 173  
www.andratec.com