

cNovate

M E D I C A L

Catalogue

- ▀ Peripheral Balloon Dilatation Catheter
- ▀ Coronary Balloon Dilatation Catheter

Thrift

Semi-Compliant 0.035 OTW PTA

2

Millau NC

Non-Compliant 0.035 OTW PTA

6

Skinny

Semi-Compliant 0.018 OTW PTA

8

Erasmus NC

Non-Compliant 0.018 OTW PTA

14

Zeeland

Semi-Compliant 0.014 OTW PTA

18

Rialto NC

Non-Compliant 0.014 OTW PTA

22

Helix

Semi-Compliant PTCA

26

Vecchio NC

Non-Compliant PTCA

28

Thrift

PTA Balloon Catheter Semi-Compliant 0.035 OTW

- **Nominal Pressure:** 6atm
- **Rated Burst Pressure (RBP):** 14atm
- **Average Burst Pressure:** 20atm
- **Sheath Compatibility:** 5F (Ø3.0-5.0mm); 6F (Ø6.0-8.0mm); 7F (Ø9.0-10.0mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
Catheter Working Length: 40cm				819-060020-040	20mm	6.0mm	40cm	819-090040-040	40mm	9.0mm	40cm
819-030020-040	20mm	3.0mm	40cm	819-060030-040	30mm	6.0mm	40cm	819-090060-040	60mm	9.0mm	40cm
819-030030-040	30mm	3.0mm	40cm	819-060040-040	40mm	6.0mm	40cm	819-090080-040	80mm	9.0mm	40cm
819-030040-040	40mm	3.0mm	40cm	819-060060-040	60mm	6.0mm	40cm	819-090100-040	100mm	9.0mm	40cm
819-030060-040	60mm	3.0mm	40cm	819-060080-040	80mm	6.0mm	40cm	819-090120-040	120mm	9.0mm	40cm
819-030080-040	80mm	3.0mm	40cm	819-060100-040	100mm	6.0mm	40cm	819-100020-040	20mm	10.0mm	40cm
819-030100-040	100mm	3.0mm	40cm	819-060120-040	120mm	6.0mm	40cm	819-100030-040	30mm	10.0mm	40cm
819-030120-040	120mm	3.0mm	40cm	819-060150-040	150mm	6.0mm	40cm	819-100040-040	40mm	10.0mm	40cm
819-030150-040	150mm	3.0mm	40cm	819-060200-040	200mm	6.0mm	40cm	819-100060-040	60mm	10.0mm	40cm
819-030200-040	200mm	3.0mm	40cm	819-070020-040	20mm	7.0mm	40cm	819-100080-040	80mm	10.0mm	40cm
819-040020-040	20mm	4.0mm	40cm	819-070030-040	30mm	7.0mm	40cm	Catheter Working Length: 70cm			
819-040030-040	30mm	4.0mm	40cm	819-070040-040	40mm	7.0mm	40cm	819-030020-070	20mm	3.0mm	70cm
819-040040-040	40mm	4.0mm	40cm	819-070060-040	60mm	7.0mm	40cm	819-030030-070	30mm	3.0mm	70cm
819-040060-040	60mm	4.0mm	40cm	819-070080-040	80mm	7.0mm	40cm	819-030040-070	40mm	3.0mm	70cm
819-040080-040	80mm	4.0mm	40cm	819-070100-040	100mm	7.0mm	40cm	819-030060-070	60mm	3.0mm	70cm
819-040100-040	100mm	4.0mm	40cm	819-070120-040	120mm	7.0mm	40cm	819-030080-070	80mm	3.0mm	70cm
819-040120-040	120mm	4.0mm	40cm	819-070150-040	150mm	7.0mm	40cm	819-030100-070	100mm	3.0mm	70cm
819-040150-040	150mm	4.0mm	40cm	819-070200-040	200mm	7.0mm	40cm	819-030120-070	120mm	3.0mm	70cm
819-040200-040	200mm	4.0mm	40cm	819-080020-040	20mm	8.0mm	40cm	819-030150-070	150mm	3.0mm	70cm
819-050020-040	20mm	5.0mm	40cm	819-080030-040	30mm	8.0mm	40cm	819-030200-070	200mm	3.0mm	70cm
819-050030-040	30mm	5.0mm	40cm	819-080040-040	40mm	8.0mm	40cm	819-040020-070	20mm	4.0mm	70cm
819-050040-040	40mm	5.0mm	40cm	819-080060-040	60mm	8.0mm	40cm	819-040030-070	30mm	4.0mm	70cm
819-050060-040	60mm	5.0mm	40cm	819-080080-040	80mm	8.0mm	40cm	819-040040-070	40mm	4.0mm	70cm
819-050080-040	80mm	5.0mm	40cm	819-080100-040	100mm	8.0mm	40cm	819-040060-070	60mm	4.0mm	70cm
819-050100-040	100mm	5.0mm	40cm	819-080120-040	120mm	8.0mm	40cm	819-040080-070	80mm	4.0mm	70cm
819-050120-040	120mm	5.0mm	40cm	819-080150-040	150mm	8.0mm	40cm	819-040100-070	100mm	4.0mm	70cm
819-050150-040	150mm	5.0mm	40cm	819-090020-040	20mm	9.0mm	40cm	819-040120-070	120mm	4.0mm	70cm
819-050200-040	200mm	5.0mm	40cm	819-090030-040	30mm	9.0mm	40cm				

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
819-040150-070	150mm	4.0mm	70cm	819-100040-070	40mm	10.0mm	70cm	819-070100-090	100mm	7.0mm	90cm
819-040200-070	200mm	4.0mm	70cm	819-100060-070	60mm	10.0mm	70cm	819-070120-090	120mm	7.0mm	90cm
819-050020-070	20mm	5.0mm	70cm	819-100080-070	80mm	10.0mm	70cm	819-070150-090	150mm	7.0mm	90cm
819-050030-070	30mm	5.0mm	70cm	Catheter Working Length: 90cm			819-070200-090	200mm	7.0mm	90cm	
819-050040-070	40mm	5.0mm	70cm	819-030020-090	20mm	3.0mm	90cm	819-080020-090	20mm	8.0mm	90cm
819-050060-070	60mm	5.0mm	70cm	819-030030-090	30mm	3.0mm	90cm	819-080030-090	30mm	8.0mm	90cm
819-050080-070	80mm	5.0mm	70cm	819-030040-090	40mm	3.0mm	90cm	819-080040-090	40mm	8.0mm	90cm
819-050100-070	100mm	5.0mm	70cm	819-030060-090	60mm	3.0mm	90cm	819-080060-090	60mm	8.0mm	90cm
819-050120-070	120mm	5.0mm	70cm	819-030080-090	80mm	3.0mm	90cm	819-080080-090	80mm	8.0mm	90cm
819-050150-070	150mm	5.0mm	70cm	819-030100-090	100mm	3.0mm	90cm	819-080100-090	100mm	8.0mm	90cm
819-050200-070	200mm	5.0mm	70cm	819-030120-090	120mm	3.0mm	90cm	819-080120-090	120mm	8.0mm	90cm
819-060020-070	20mm	6.0mm	70cm	819-030150-090	150mm	3.0mm	90cm	819-080150-090	150mm	8.0mm	90cm
819-060030-070	30mm	6.0mm	70cm	819-030200-090	200mm	3.0mm	90cm	819-090020-090	20mm	9.0mm	90cm
819-060040-070	40mm	6.0mm	70cm	819-040020-090	20mm	4.0mm	90cm	819-090030-090	30mm	9.0mm	90cm
819-060060-070	60mm	6.0mm	70cm	819-040030-090	30mm	4.0mm	90cm	819-090040-090	40mm	9.0mm	90cm
819-060080-070	80mm	6.0mm	70cm	819-040040-090	40mm	4.0mm	90cm	819-090060-090	60mm	9.0mm	90cm
819-060100-070	100mm	6.0mm	70cm	819-040060-090	60mm	4.0mm	90cm	819-090080-090	80mm	9.0mm	90cm
819-060120-070	120mm	6.0mm	70cm	819-040080-090	80mm	4.0mm	90cm	819-090100-090	100mm	9.0mm	90cm
819-060150-070	150mm	6.0mm	70cm	819-040100-090	100mm	4.0mm	90cm	819-090120-090	120mm	9.0mm	90cm
819-060200-070	200mm	6.0mm	70cm	819-040120-090	120mm	4.0mm	90cm	819-100020-090	20mm	10.0mm	90cm
819-070020-070	20mm	7.0mm	70cm	819-040150-090	150mm	4.0mm	90cm	819-100030-090	30mm	10.0mm	90cm
819-070030-070	30mm	7.0mm	70cm	819-040200-090	200mm	4.0mm	90cm	819-100040-090	40mm	10.0mm	90cm
819-070040-070	40mm	7.0mm	70cm	819-050020-090	20mm	5.0mm	90cm	819-100060-090	60mm	10.0mm	90cm
819-070060-070	60mm	7.0mm	70cm	819-050030-090	30mm	5.0mm	90cm	819-100080-090	80mm	10.0mm	90cm
819-070080-070	80mm	7.0mm	70cm	819-050040-090	40mm	5.0mm	90cm	Catheter Working Length: 150cm			
819-070100-070	100mm	7.0mm	70cm	819-050060-090	60mm	5.0mm	90cm	819-030020-150	20mm	3.0mm	150cm
819-070120-070	120mm	7.0mm	70cm	819-050080-090	80mm	5.0mm	90cm	819-030030-150	30mm	3.0mm	150cm
819-070150-070	150mm	7.0mm	70cm	819-050100-090	100mm	5.0mm	90cm	819-030040-150	40mm	3.0mm	150cm
819-070200-070	200mm	7.0mm	70cm	819-050120-090	120mm	5.0mm	90cm	819-030060-150	60mm	3.0mm	150cm
819-080020-070	20mm	8.0mm	70cm	819-050150-090	150mm	5.0mm	90cm	819-030080-150	80mm	3.0mm	150cm
819-080030-070	30mm	8.0mm	70cm	819-050200-090	200mm	5.0mm	90cm	819-030100-150	100mm	3.0mm	150cm
819-080040-070	40mm	8.0mm	70cm	819-060020-090	20mm	6.0mm	90cm	819-030120-150	120mm	3.0mm	150cm
819-080060-070	60mm	8.0mm	70cm	819-060030-090	30mm	6.0mm	90cm	819-030150-150	150mm	3.0mm	150cm
819-080080-070	80mm	8.0mm	70cm	819-060040-090	40mm	6.0mm	90cm	819-030200-150	200mm	3.0mm	150cm
819-080100-070	100mm	8.0mm	70cm	819-060060-090	60mm	6.0mm	90cm	819-040020-150	20mm	4.0mm	150cm
819-080120-070	120mm	8.0mm	70cm	819-060080-090	80mm	6.0mm	90cm	819-040030-150	30mm	4.0mm	150cm
819-080150-070	150mm	8.0mm	70cm	819-060100-090	100mm	6.0mm	90cm	819-040040-150	40mm	4.0mm	150cm
819-090020-070	20mm	9.0mm	70cm	819-060120-090	120mm	6.0mm	90cm	819-040060-150	60mm	4.0mm	150cm
819-090030-070	30mm	9.0mm	70cm	819-060150-090	150mm	6.0mm	90cm	819-040080-150	80mm	4.0mm	150cm
819-090040-070	40mm	9.0mm	70cm	819-060200-090	200mm	6.0mm	90cm	819-040100-150	100mm	4.0mm	150cm
819-090060-070	60mm	9.0mm	70cm	819-070020-090	20mm	7.0mm	90cm	819-040120-150	120mm	4.0mm	150cm
819-090080-070	80mm	9.0mm	70cm	819-070030-090	30mm	7.0mm	90cm	819-040150-150	150mm	4.0mm	150cm
819-090100-070	100mm	9.0mm	70cm	819-070040-090	40mm	7.0mm	90cm	819-040200-150	200mm	4.0mm	150cm
819-090120-070	120mm	9.0mm	70cm	819-070060-090	60mm	7.0mm	90cm	819-050020-150	20mm	5.0mm	150cm
819-100020-070	20mm	10.0mm	70cm	819-070080-090	80mm	7.0mm	90cm	819-050030-150	30mm	5.0mm	150cm
819-100030-070	30mm	10.0mm	70cm								

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
819-050040-150	40mm	5.0mm	150cm	819-070020-150	20mm	7.0mm	150cm	819-080150-150	150mm	8.0mm	150cm
819-050060-150	60mm	5.0mm	150cm	819-070030-150	30mm	7.0mm	150cm	819-090020-150	20mm	9.0mm	150cm
819-050080-150	80mm	5.0mm	150cm	819-070040-150	40mm	7.0mm	150cm	819-090030-150	30mm	9.0mm	150cm
819-050100-150	100mm	5.0mm	150cm	819-070060-150	60mm	7.0mm	150cm	819-090040-150	40mm	9.0mm	150cm
819-050120-150	120mm	5.0mm	150cm	819-070080-150	80mm	7.0mm	150cm	819-090060-150	60mm	9.0mm	150cm
819-050150-150	150mm	5.0mm	150cm	819-070100-150	100mm	7.0mm	150cm	819-090080-150	80mm	9.0mm	150cm
819-050200-150	200mm	5.0mm	150cm	819-070120-150	120mm	7.0mm	150cm	819-090100-150	100mm	9.0mm	150cm
819-060020-150	20mm	6.0mm	150cm	819-070150-150	150mm	7.0mm	150cm	819-090120-150	120mm	9.0mm	150cm
819-060030-150	30mm	6.0mm	150cm	819-070200-150	200mm	7.0mm	150cm	819-100020-150	20mm	10.0mm	150cm
819-060040-150	40mm	6.0mm	150cm	819-080020-150	20mm	8.0mm	150cm	819-100030-150	30mm	10.0mm	150cm
819-060060-150	60mm	6.0mm	150cm	819-080030-150	30mm	8.0mm	150cm	819-100040-150	40mm	10.0mm	150cm
819-060080-150	80mm	6.0mm	150cm	819-080040-150	40mm	8.0mm	150cm	819-100060-150	60mm	10.0mm	150cm
819-060100-150	100mm	6.0mm	150cm	819-080060-150	60mm	8.0mm	150cm	819-100080-150	80mm	10.0mm	150cm
819-060120-150	120mm	6.0mm	150cm	819-080080-150	80mm	8.0mm	150cm				
819-060150-150	150mm	6.0mm	150cm	819-080100-150	100mm	8.0mm	150cm				
819-060200-150	200mm	6.0mm	150cm	819-080120-150	120mm	8.0mm	150cm				





Millau NC

PTA Balloon Catheter Non-Compliant 0.035 OTW

- **Nominal Pressure:** 12atm
- **Rated Burst Pressure (RBP):** 22atm (Ø3.0-4.0mm); 20atm (Ø5.0-6.0mm); 18atm (Ø7.0-10.0mm)
- **Average Burst Pressure:** 30atm (Ø3.0-5.0mm); 26atm (Ø6.0-10.0mm)
- **Sheath Compatibility:** 5F (Ø3.0-5.0mm); 6F (Ø6.0-8.0mm); 7F (Ø9.0-10.0mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
Catheter Working Length: 40cm				Catheter Working Length: 70cm							
820-030020-040	20mm	3.0mm	40cm	820-060150-040	150mm	6.0mm	40cm	820-030020-070	20mm	3.0mm	70cm
820-030030-040	30mm	3.0mm	40cm	820-070020-040	20mm	7.0mm	40cm	820-030030-070	30mm	3.0mm	70cm
820-030040-040	40mm	3.0mm	40cm	820-070030-040	30mm	7.0mm	40cm	820-030040-070	40mm	3.0mm	70cm
820-030060-040	60mm	3.0mm	40cm	820-070040-040	40mm	7.0mm	40cm	820-030060-070	60mm	3.0mm	70cm
820-030080-040	80mm	3.0mm	40cm	820-070060-040	60mm	7.0mm	40cm	820-030080-070	80mm	3.0mm	70cm
820-030100-040	100mm	3.0mm	40cm	820-070080-040	80mm	7.0mm	40cm	820-030100-070	100mm	3.0mm	70cm
820-030150-040	150mm	3.0mm	40cm	820-070100-040	100mm	7.0mm	40cm	820-030150-070	150mm	3.0mm	70cm
820-040020-040	20mm	4.0mm	40cm	820-070150-040	150mm	7.0mm	40cm	820-040020-070	20mm	4.0mm	70cm
820-040030-040	30mm	4.0mm	40cm	820-080020-040	20mm	8.0mm	40cm	820-040030-070	30mm	4.0mm	70cm
820-040040-040	40mm	4.0mm	40cm	820-080030-040	30mm	8.0mm	40cm	820-040040-070	40mm	4.0mm	70cm
820-040060-040	60mm	4.0mm	40cm	820-080040-040	40mm	8.0mm	40cm	820-040060-070	60mm	4.0mm	70cm
820-040080-040	80mm	4.0mm	40cm	820-080060-040	60mm	8.0mm	40cm	820-040080-070	80mm	4.0mm	70cm
820-040100-040	100mm	4.0mm	40cm	820-080080-040	80mm	8.0mm	40cm	820-040100-070	100mm	4.0mm	70cm
820-040150-040	150mm	4.0mm	40cm	820-080100-040	100mm	8.0mm	40cm	820-040150-070	150mm	4.0mm	70cm
820-050020-040	20mm	5.0mm	40cm	820-080150-040	150mm	8.0mm	40cm	820-050020-070	20mm	5.0mm	70cm
820-050030-040	30mm	5.0mm	40cm	820-090020-040	20mm	9.0mm	40cm	820-050030-070	30mm	5.0mm	70cm
820-050040-040	40mm	5.0mm	40cm	820-090030-040	30mm	9.0mm	40cm	820-050040-070	40mm	5.0mm	70cm
820-050060-040	60mm	5.0mm	40cm	820-090040-040	40mm	9.0mm	40cm	820-050060-070	60mm	5.0mm	70cm
820-050080-040	80mm	5.0mm	40cm	820-090060-040	60mm	9.0mm	40cm	820-050080-070	80mm	5.0mm	70cm
820-050100-040	100mm	5.0mm	40cm	820-090080-040	80mm	9.0mm	40cm	820-050100-070	100mm	5.0mm	70cm
820-050150-040	150mm	5.0mm	40cm	820-090100-040	100mm	9.0mm	40cm	820-050150-070	150mm	5.0mm	70cm
820-060020-040	20mm	6.0mm	40cm	820-100020-040	20mm	10.0mm	40cm	820-060020-070	20mm	6.0mm	70cm
820-060030-040	30mm	6.0mm	40cm	820-100030-040	30mm	10.0mm	40cm	820-060030-070	30mm	6.0mm	70cm
820-060040-040	40mm	6.0mm	40cm	820-100040-040	40mm	10.0mm	40cm	820-060040-070	40mm	6.0mm	70cm
820-060060-040	60mm	6.0mm	40cm	820-100060-040	60mm	10.0mm	40cm	820-060060-070	60mm	6.0mm	70cm
820-060080-040	80mm	6.0mm	40cm	820-100080-040	80mm	10.0mm	40cm	820-060080-070	80mm	6.0mm	70cm
820-060100-040	100mm	6.0mm	40cm					820-060100-070	100mm	6.0mm	70cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
820-060150-070	150mm	6.0mm	70cm	820-050080-090	80mm	5.0mm	90cm	820-040040-150	40mm	4.0mm	150cm
820-070020-070	20mm	7.0mm	70cm	820-050100-090	100mm	5.0mm	90cm	820-040060-150	60mm	4.0mm	150cm
820-070030-070	30mm	7.0mm	70cm	820-050150-090	150mm	5.0mm	90cm	820-040080-150	80mm	4.0mm	150cm
820-070040-070	40mm	7.0mm	70cm	820-060020-090	20mm	6.0mm	90cm	820-040100-150	100mm	4.0mm	150cm
820-070060-070	60mm	7.0mm	70cm	820-060030-090	30mm	6.0mm	90cm	820-040150-150	150mm	4.0mm	150cm
820-070080-070	80mm	7.0mm	70cm	820-060040-090	40mm	6.0mm	90cm	820-050020-150	20mm	5.0mm	150cm
820-070100-070	100mm	7.0mm	70cm	820-060060-090	60mm	6.0mm	90cm	820-050030-150	30mm	5.0mm	150cm
820-070150-070	150mm	7.0mm	70cm	820-060080-090	80mm	6.0mm	90cm	820-050040-150	40mm	5.0mm	150cm
820-080020-070	20mm	8.0mm	70cm	820-060100-090	100mm	6.0mm	90cm	820-050060-150	60mm	5.0mm	150cm
820-080030-070	30mm	8.0mm	70cm	820-060150-090	150mm	6.0mm	90cm	820-050080-150	80mm	5.0mm	150cm
820-080040-070	40mm	8.0mm	70cm	820-070020-090	20mm	7.0mm	90cm	820-050100-150	100mm	5.0mm	150cm
820-080060-070	60mm	8.0mm	70cm	820-070030-090	30mm	7.0mm	90cm	820-050150-150	150mm	5.0mm	150cm
820-080080-070	80mm	8.0mm	70cm	820-070040-090	40mm	7.0mm	90cm	820-060020-150	20mm	6.0mm	150cm
820-080100-070	100mm	8.0mm	70cm	820-070060-090	60mm	7.0mm	90cm	820-060030-150	30mm	6.0mm	150cm
820-080150-070	150mm	8.0mm	70cm	820-070080-090	80mm	7.0mm	90cm	820-060040-150	40mm	6.0mm	150cm
820-090020-070	20mm	9.0mm	70cm	820-070100-090	100mm	7.0mm	90cm	820-060060-150	60mm	6.0mm	150cm
820-090030-070	30mm	9.0mm	70cm	820-070150-090	150mm	7.0mm	90cm	820-060080-150	80mm	6.0mm	150cm
820-090040-070	40mm	9.0mm	70cm	820-080020-090	20mm	8.0mm	90cm	820-060100-150	100mm	6.0mm	150cm
820-090060-070	60mm	9.0mm	70cm	820-080030-090	30mm	8.0mm	90cm	820-060150-150	150mm	6.0mm	150cm
820-090080-070	80mm	9.0mm	70cm	820-080040-090	40mm	8.0mm	90cm	820-070020-150	20mm	7.0mm	150cm
820-090100-070	100mm	9.0mm	70cm	820-080060-090	60mm	8.0mm	90cm	820-070030-150	30mm	7.0mm	150cm
820-100020-070	20mm	10.0mm	70cm	820-080080-090	80mm	8.0mm	90cm	820-070040-150	40mm	7.0mm	150cm
820-100030-070	30mm	10.0mm	70cm	820-080100-090	100mm	8.0mm	90cm	820-070060-150	60mm	7.0mm	150cm
820-100040-070	40mm	10.0mm	70cm	820-080150-090	150mm	8.0mm	90cm	820-070080-150	80mm	7.0mm	150cm
820-100060-070	60mm	10.0mm	70cm	820-090020-090	20mm	9.0mm	90cm	820-070100-150	100mm	7.0mm	150cm
820-100080-070	80mm	10.0mm	70cm	820-090030-090	30mm	9.0mm	90cm	820-070150-150	150mm	7.0mm	150cm
Catheter Working Length: 90cm				820-090040-090	40mm	9.0mm	90cm	820-080020-150	20mm	8.0mm	150cm
820-030020-090	20mm	3.0mm	90cm	820-090060-090	60mm	9.0mm	90cm	820-080030-150	30mm	8.0mm	150cm
820-030030-090	30mm	3.0mm	90cm	820-090080-090	80mm	9.0mm	90cm	820-080040-150	40mm	8.0mm	150cm
820-030040-090	40mm	3.0mm	90cm	820-090100-090	100mm	9.0mm	90cm	820-080060-150	60mm	8.0mm	150cm
820-030060-090	60mm	3.0mm	90cm	820-100020-090	20mm	10.0mm	90cm	820-080080-150	80mm	8.0mm	150cm
820-030080-090	80mm	3.0mm	90cm	820-100030-090	30mm	10.0mm	90cm	820-080100-150	100mm	8.0mm	150cm
820-030100-090	100mm	3.0mm	90cm	820-100040-090	40mm	10.0mm	90cm	820-080150-150	150mm	8.0mm	150cm
820-030150-090	150mm	3.0mm	90cm	820-100060-090	60mm	10.0mm	90cm	820-090020-150	20mm	9.0mm	150cm
820-040020-090	20mm	4.0mm	90cm	820-100080-090	80mm	10.0mm	90cm	820-090030-150	30mm	9.0mm	150cm
820-040030-090	30mm	4.0mm	90cm	Catheter Working Length: 150cm				820-090040-150	40mm	9.0mm	150cm
820-040040-090	40mm	4.0mm	90cm	820-030020-150	20mm	3.0mm	150cm	820-090060-150	60mm	9.0mm	150cm
820-040060-090	60mm	4.0mm	90cm	820-030030-150	30mm	3.0mm	150cm	820-090080-150	80mm	9.0mm	150cm
820-040080-090	80mm	4.0mm	90cm	820-030040-150	40mm	3.0mm	150cm	820-090100-150	100mm	9.0mm	150cm
820-040100-090	100mm	4.0mm	90cm	820-030060-150	60mm	3.0mm	150cm	820-100020-150	20mm	10.0mm	150cm
820-040150-090	150mm	4.0mm	90cm	820-030080-150	80mm	3.0mm	150cm	820-100030-150	30mm	10.0mm	150cm
820-050020-090	20mm	5.0mm	90cm	820-030100-150	100mm	3.0mm	150cm	820-100040-150	40mm	10.0mm	150cm
820-050030-090	30mm	5.0mm	90cm	820-030150-150	150mm	3.0mm	150cm	820-100060-150	60mm	10.0mm	150cm
820-050040-090	40mm	5.0mm	90cm	820-040020-150	20mm	4.0mm	150cm	820-100080-150	80mm	10.0mm	150cm
820-050060-090	60mm	5.0mm	90cm	820-040030-150	30mm	4.0mm	150cm				



Skinny

PTA Balloon Catheter Semi-compliant 0.018 OTW

- **Nominal Pressure:** 6atm
- **Rated Burst Pressure (RBP):** 14atm
- **Average Burst Pressure:** 20atm
- **Sheath Compatibility:** 4F (Ø1.5-4.0mm x 5-200mm; Ø4.5-6.0 x 5-60mm)
5F (Ø4.5-6.0mm x 80-200mm; Ø7.0 x 5-200mm)
6F (Ø8.0-10.0mm x 5-200mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
Catheter Working Length: 40cm				803-025020-040	20mm	2.5mm	40cm	803-035080-040	80mm	3.5mm	40cm
803-015005-040	5mm	1.5mm	40cm	803-025030-040	30mm	2.5mm	40cm	803-035100-040	100mm	3.5mm	40cm
803-015010-040	10mm	1.5mm	40cm	803-025040-040	40mm	2.5mm	40cm	803-035120-040	120mm	3.5mm	40cm
803-015015-040	15mm	1.5mm	40cm	803-025060-040	60mm	2.5mm	40cm	803-035150-040	150mm	3.5mm	40cm
803-015020-040	20mm	1.5mm	40cm	803-025080-040	80mm	2.5mm	40cm	803-035200-040	200mm	3.5mm	40cm
803-015030-040	30mm	1.5mm	40cm	803-025100-040	100mm	2.5mm	40cm	803-040005-040	5mm	4.0mm	40cm
803-015040-040	40mm	1.5mm	40cm	803-025120-040	120mm	2.5mm	40cm	803-040010-040	10mm	4.0mm	40cm
803-015060-040	60mm	1.5mm	40cm	803-025150-040	150mm	2.5mm	40cm	803-040015-040	15mm	4.0mm	40cm
803-015080-040	80mm	1.5mm	40cm	803-025200-040	200mm	2.5mm	40cm	803-040020-040	20mm	4.0mm	40cm
803-015100-040	100mm	1.5mm	40cm	803-030005-040	5mm	3.0mm	40cm	803-040030-040	30mm	4.0mm	40cm
803-015120-040	120mm	1.5mm	40cm	803-030010-040	10mm	3.0mm	40cm	803-040040-040	40mm	4.0mm	40cm
803-015150-040	150mm	1.5mm	40cm	803-030015-040	15mm	3.0mm	40cm	803-040060-040	60mm	4.0mm	40cm
803-015200-040	200mm	1.5mm	40cm	803-030020-040	20mm	3.0mm	40cm	803-040080-040	80mm	4.0mm	40cm
803-020005-040	5mm	2.0mm	40cm	803-030030-040	30mm	3.0mm	40cm	803-040100-040	100mm	4.0mm	40cm
803-020010-040	10mm	2.0mm	40cm	803-030040-040	40mm	3.0mm	40cm	803-040120-040	120mm	4.0mm	40cm
803-020015-040	15mm	2.0mm	40cm	803-030060-040	60mm	3.0mm	40cm	803-040150-040	150mm	4.0mm	40cm
803-020020-040	20mm	2.0mm	40cm	803-030080-040	80mm	3.0mm	40cm	803-040200-040	200mm	4.0mm	40cm
803-020030-040	30mm	2.0mm	40cm	803-030100-040	100mm	3.0mm	40cm	803-045005-040	5mm	4.5mm	40cm
803-020040-040	40mm	2.0mm	40cm	803-030120-040	120mm	3.0mm	40cm	803-045010-040	10mm	4.5mm	40cm
803-020060-040	60mm	2.0mm	40cm	803-030150-040	150mm	3.0mm	40cm	803-045015-040	15mm	4.5mm	40cm
803-020080-040	80mm	2.0mm	40cm	803-030200-040	200mm	3.0mm	40cm	803-045020-040	20mm	4.5mm	40cm
803-020100-040	100mm	2.0mm	40cm	803-035005-040	5mm	3.5mm	40cm	803-045030-040	30mm	4.5mm	40cm
803-020120-040	120mm	2.0mm	40cm	803-035010-040	10mm	3.5mm	40cm	803-045040-040	40mm	4.5mm	40cm
803-020150-040	150mm	2.0mm	40cm	803-035015-040	15mm	3.5mm	40cm	803-045060-040	60mm	4.5mm	40cm
803-020200-040	200mm	2.0mm	40cm	803-035020-040	20mm	3.5mm	40cm	803-045080-040	80mm	4.5mm	40cm
803-025005-040	5mm	2.5mm	40cm	803-035030-040	30mm	3.5mm	40cm	803-045100-040	100mm	4.5mm	40cm
803-025010-040	10mm	2.5mm	40cm	803-035040-040	40mm	3.5mm	40cm	803-045120-040	120mm	4.5mm	40cm
803-025015-040	15mm	2.5mm	40cm	803-035060-040	60mm	3.5mm	40cm	803-045150-040	150mm	4.5mm	40cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
803-035150-070	150mm	3.5mm	70cm	803-055100-070	100mm	5.5mm	70cm	803-090080-070	80mm	9.0mm	70cm
803-035200-070	200mm	3.5mm	70cm	803-055120-070	120mm	5.5mm	70cm	803-090100-070	100mm	9.0mm	70cm
803-040005-070	5mm	4.0mm	70cm	803-055150-070	150mm	5.5mm	70cm	803-090120-070	120mm	9.0mm	70cm
803-040010-070	10mm	4.0mm	70cm	803-055200-070	200mm	5.5mm	70cm	803-100005-070	5mm	8.0mm	70cm
803-040015-070	15mm	4.0mm	70cm	803-060005-070	5mm	6.0mm	70cm	803-100010-070	10mm	10.0mm	70cm
803-040020-070	20mm	4.0mm	70cm	803-060010-070	10mm	6.0mm	70cm	803-100015-070	15mm	10.0mm	70cm
803-040030-070	30mm	4.0mm	70cm	803-060015-070	15mm	6.0mm	70cm	803-100020-070	20mm	10.0mm	70cm
803-040040-070	40mm	4.0mm	70cm	803-060020-070	20mm	6.0mm	70cm	803-100030-070	30mm	10.0mm	70cm
803-040060-070	60mm	4.0mm	70cm	803-060030-070	30mm	6.0mm	70cm	803-100040-070	40mm	10.0mm	70cm
803-040080-070	80mm	4.0mm	70cm	803-060040-070	40mm	6.0mm	70cm	803-100060-070	60mm	10.0mm	70cm
803-040100-070	100mm	4.0mm	70cm	803-060060-070	60mm	6.0mm	70cm	803-100080-070	80mm	10.0mm	70cm
803-040120-070	120mm	4.0mm	70cm	803-060080-070	80mm	6.0mm	70cm				
803-040150-070	150mm	4.0mm	70cm	803-060100-070	100mm	6.0mm	70cm	Catheter Working Length: 90cm			
803-040200-070	200mm	4.0mm	70cm	803-060120-070	120mm	6.0mm	70cm	803-015005-090	5mm	1.5mm	90cm
803-045005-070	5mm	4.5mm	70cm	803-060150-070	150mm	6.0mm	70cm	803-015010-090	10mm	1.5mm	90cm
803-045010-070	10mm	4.5mm	70cm	803-060200-070	200mm	6.0mm	70cm	803-015015-090	15mm	1.5mm	90cm
803-045015-070	15mm	4.5mm	70cm	803-070005-070	5mm	7.0mm	70cm	803-015020-090	20mm	1.5mm	90cm
803-045020-070	20mm	4.5mm	70cm	803-070010-070	10mm	7.0mm	70cm	803-015030-090	30mm	1.5mm	90cm
803-045030-070	30mm	4.5mm	70cm	803-070015-070	15mm	7.0mm	70cm	803-015040-090	40mm	1.5mm	90cm
803-045040-070	40mm	4.5mm	70cm	803-070020-070	20mm	7.0mm	70cm	803-015060-090	60mm	1.5mm	90cm
803-045060-070	60mm	4.5mm	70cm	803-070030-070	30mm	7.0mm	70cm	803-015080-090	80mm	1.5mm	90cm
803-045080-070	80mm	4.5mm	70cm	803-070040-070	40mm	7.0mm	70cm	803-015100-090	100mm	1.5mm	90cm
803-045100-070	100mm	4.5mm	70cm	803-070060-070	60mm	7.0mm	70cm	803-015120-090	120mm	1.5mm	90cm
803-045120-070	120mm	4.5mm	70cm	803-070080-070	80mm	7.0mm	70cm	803-015150-090	150mm	1.5mm	90cm
803-045150-070	150mm	4.5mm	70cm	803-070100-070	100mm	7.0mm	70cm	803-015200-090	200mm	1.5mm	90cm
803-045200-070	200mm	4.5mm	70cm	803-070120-070	120mm	7.0mm	70cm	803-020005-090	5mm	2.0mm	90cm
803-050005-070	5mm	5.0mm	70cm	803-070150-070	150mm	7.0mm	70cm	803-020010-090	10mm	2.0mm	90cm
803-050010-070	10mm	5.0mm	70cm	803-070200-070	200mm	7.0mm	70cm	803-020015-090	15mm	2.0mm	90cm
803-050015-070	15mm	5.0mm	70cm	803-080005-070	5mm	8.0mm	70cm	803-020020-090	20mm	2.0mm	90cm
803-050020-070	20mm	5.0mm	70cm	803-080010-070	10mm	8.0mm	70cm	803-020030-090	30mm	2.0mm	90cm
803-050030-070	30mm	5.0mm	70cm	803-080015-070	15mm	8.0mm	70cm	803-020040-090	40mm	2.0mm	90cm
803-050040-070	40mm	5.0mm	70cm	803-080020-070	20mm	8.0mm	70cm	803-020060-090	60mm	2.0mm	90cm
803-050060-070	60mm	5.0mm	70cm	803-080030-070	30mm	8.0mm	70cm	803-020080-090	80mm	2.0mm	90cm
803-050080-070	80mm	5.0mm	70cm	803-080040-070	40mm	8.0mm	70cm	803-020100-090	100mm	2.0mm	90cm
803-050100-070	100mm	5.0mm	70cm	803-080060-070	60mm	8.0mm	70cm	803-020120-090	120mm	2.0mm	90cm
803-050120-070	120mm	5.0mm	70cm	803-080080-070	80mm	8.0mm	70cm	803-020150-090	150mm	2.0mm	90cm
803-050150-070	150mm	5.0mm	70cm	803-080100-070	100mm	8.0mm	70cm	803-020200-090	200mm	2.0mm	90cm
803-050200-070	200mm	5.0mm	70cm	803-080120-070	120mm	8.0mm	70cm	803-025005-090	5mm	2.5mm	90cm
803-055005-070	5mm	5.5mm	70cm	803-080150-070	150mm	8.0mm	70cm	803-025010-090	10mm	2.5mm	90cm
803-055010-070	10mm	5.5mm	70cm	803-090005-070	5mm	9.0mm	70cm	803-025015-090	15mm	2.5mm	90cm
803-055015-070	15mm	5.5mm	70cm	803-090010-070	10mm	9.0mm	70cm	803-025020-090	20mm	2.5mm	90cm
803-055020-070	20mm	5.5mm	70cm	803-090015-070	15mm	9.0mm	70cm	803-025030-090	30mm	2.5mm	90cm
803-055030-070	30mm	5.5mm	70cm	803-090020-070	20mm	9.0mm	70cm	803-025040-090	40mm	2.5mm	90cm
803-055040-070	40mm	5.5mm	70cm	803-090030-070	30mm	9.0mm	70cm	803-025060-090	60mm	2.5mm	90cm
803-055060-070	60mm	5.5mm	70cm	803-090040-070	40mm	9.0mm	70cm	803-025080-090	80mm	2.5mm	90cm
803-055080-070	80mm	5.5mm	70cm	803-090060-070	60mm	9.0mm	70cm	803-025100-090	100mm	2.5mm	90cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
803-015100-150	100mm	1.5mm	150cm	803-035060-150	60mm	3.5mm	150cm	803-055030-150	30mm	5.5mm	150cm
803-015120-150	120mm	1.5mm	150cm	803-035080-150	80mm	3.5mm	150cm	803-055040-150	40mm	5.5mm	150cm
803-015150-150	150mm	1.5mm	150cm	803-035100-150	100mm	3.5mm	150cm	803-055060-150	60mm	5.5mm	150cm
803-015200-150	200mm	1.5mm	150cm	803-035120-150	120mm	3.5mm	150cm	803-055080-150	80mm	5.5mm	150cm
803-020005-150	5mm	2.0mm	150cm	803-035150-150	150mm	3.5mm	150cm	803-055100-150	100mm	5.5mm	150cm
803-020010-150	10mm	2.0mm	150cm	803-035200-150	200mm	3.5mm	150cm	803-055120-150	120mm	5.5mm	150cm
803-020015-150	15mm	2.0mm	150cm	803-040005-150	5mm	4.0mm	150cm	803-055150-150	150mm	5.5mm	150cm
803-020020-150	20mm	2.0mm	150cm	803-040010-150	10mm	4.0mm	150cm	803-055200-150	200mm	5.5mm	150cm
803-020030-150	30mm	2.0mm	150cm	803-040015-150	15mm	4.0mm	150cm	803-060005-150	5mm	6.0mm	150cm
803-020040-150	40mm	2.0mm	150cm	803-040020-150	20mm	4.0mm	150cm	803-060010-150	10mm	6.0mm	150cm
803-020060-150	60mm	2.0mm	150cm	803-040030-150	30mm	4.0mm	150cm	803-060015-150	15mm	6.0mm	150cm
803-020080-150	80mm	2.0mm	150cm	803-040040-150	40mm	4.0mm	150cm	803-060020-150	20mm	6.0mm	150cm
803-020100-150	100mm	2.0mm	150cm	803-040060-150	60mm	4.0mm	150cm	803-060030-150	30mm	6.0mm	150cm
803-020120-150	120mm	2.0mm	150cm	803-040080-150	80mm	4.0mm	150cm	803-060040-150	40mm	6.0mm	150cm
803-020150-150	150mm	2.0mm	150cm	803-040100-150	100mm	4.0mm	150cm	803-060060-150	60mm	6.0mm	150cm
803-020200-150	200mm	2.0mm	150cm	803-040120-150	120mm	4.0mm	150cm	803-060080-150	80mm	6.0mm	150cm
803-025005-150	5mm	2.5mm	150cm	803-040150-150	150mm	4.0mm	150cm	803-060100-150	100mm	6.0mm	150cm
803-025010-150	10mm	2.5mm	150cm	803-040200-150	200mm	4.0mm	150cm	803-060120-150	120mm	6.0mm	150cm
803-025015-150	15mm	2.5mm	150cm	803-045005-150	5mm	4.5mm	150cm	803-060150-150	150mm	6.0mm	150cm
803-025020-150	20mm	2.5mm	150cm	803-045010-150	10mm	4.5mm	150cm	803-060200-150	200mm	6.0mm	150cm
803-025030-150	30mm	2.5mm	150cm	803-045015-150	15mm	4.5mm	150cm	803-070005-150	5mm	7.0mm	150cm
803-025040-150	40mm	2.5mm	150cm	803-045020-150	20mm	4.5mm	150cm	803-070010-150	10mm	7.0mm	150cm
803-025060-150	60mm	2.5mm	150cm	803-045030-150	30mm	4.5mm	150cm	803-070015-150	15mm	7.0mm	150cm
803-025080-150	80mm	2.5mm	150cm	803-045040-150	40mm	4.5mm	150cm	803-070020-150	20mm	7.0mm	150cm
803-025100-150	100mm	2.5mm	150cm	803-045060-150	60mm	4.5mm	150cm	803-070030-150	30mm	7.0mm	150cm
803-025120-150	120mm	2.5mm	150cm	803-045080-150	80mm	4.5mm	150cm	803-070040-150	40mm	7.0mm	150cm
803-025150-150	150mm	2.5mm	150cm	803-045100-150	100mm	4.5mm	150cm	803-070060-150	60mm	7.0mm	150cm
803-025200-150	200mm	2.5mm	150cm	803-045120-150	120mm	4.5mm	150cm	803-070080-150	80mm	7.0mm	150cm
803-030005-150	5mm	3.0mm	150cm	803-045150-150	150mm	4.5mm	150cm	803-070100-150	100mm	7.0mm	150cm
803-030010-150	10mm	3.0mm	150cm	803-045200-150	200mm	4.5mm	150cm	803-070120-150	120mm	7.0mm	150cm
803-030015-150	15mm	3.0mm	150cm	803-050005-150	5mm	5.0mm	150cm	803-070150-150	150mm	7.0mm	150cm
803-030020-150	20mm	3.0mm	150cm	803-050010-150	10mm	5.0mm	150cm	803-070200-150	200mm	7.0mm	150cm
803-030030-150	30mm	3.0mm	150cm	803-050015-150	15mm	5.0mm	150cm	803-080005-150	5mm	8.0mm	150cm
803-030040-150	40mm	3.0mm	150cm	803-050020-150	20mm	5.0mm	150cm	803-080010-150	10mm	8.0mm	150cm
803-030060-150	60mm	3.0mm	150cm	803-050030-150	30mm	5.0mm	150cm	803-080015-150	15mm	8.0mm	150cm
803-030080-150	80mm	3.0mm	150cm	803-050040-150	40mm	5.0mm	150cm	803-080020-150	20mm	8.0mm	150cm
803-030100-150	100mm	3.0mm	150cm	803-050060-150	60mm	5.0mm	150cm	803-080030-150	30mm	8.0mm	150cm
803-030120-150	120mm	3.0mm	150cm	803-050080-150	80mm	5.0mm	150cm	803-080040-150	40mm	8.0mm	150cm
803-030150-150	150mm	3.0mm	150cm	803-050100-150	100mm	5.0mm	150cm	803-080060-150	60mm	8.0mm	150cm
803-030200-150	200mm	3.0mm	150cm	803-050120-150	120mm	5.0mm	150cm	803-080080-150	80mm	8.0mm	150cm
803-035005-150	5mm	3.5mm	150cm	803-050150-150	150mm	5.0mm	150cm	803-080100-150	100mm	8.0mm	150cm
803-035010-150	10mm	3.5mm	150cm	803-050200-150	200mm	5.0mm	150cm	803-080120-150	120mm	8.0mm	150cm
803-035015-150	15mm	3.5mm	150cm	803-055005-150	5mm	5.5mm	150cm	803-080150-150	150mm	8.0mm	150cm
803-035020-150	20mm	3.5mm	150cm	803-055010-150	10mm	5.5mm	150cm	803-090005-150	5mm	9.0mm	150cm
803-035030-150	30mm	3.5mm	150cm	803-055015-150	15mm	5.5mm	150cm	803-090010-150	10mm	9.0mm	150cm
803-035040-150	40mm	3.5mm	150cm	803-055020-150	20mm	5.5mm	150cm	803-090015-150	15mm	9.0mm	150cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
803-090020-150	20mm	9.0mm	150cm	803-090120-150	120mm	9.0mm	150cm	803-100040-150	40mm	10.0mm	150cm
803-090030-150	30mm	9.0mm	150cm	803-100005-150	5mm	8.0mm	150cm	803-100060-150	60mm	10.0mm	150cm
803-090040-150	40mm	9.0mm	150cm	803-100010-150	10mm	10.0mm	150cm	803-100080-150	80mm	10.0mm	150cm
803-090060-150	60mm	9.0mm	150cm	803-100015-150	15mm	10.0mm	150cm				
803-090080-150	80mm	9.0mm	150cm	803-100020-150	20mm	10.0mm	150cm				
803-090100-150	100mm	9.0mm	150cm	803-100030-150	30mm	10.0mm	150cm				





Erasmus NC

PTA Balloon Catheter Non-Compliant 0.018 OTW

- **Nominal Pressure:** 12atm
- **Rated Burst Pressure (RBP):** 22atm (Ø2.0-4.0mm); 20atm (Ø5.0-6.0mm); 18atm (Ø7.0-8.0mm)
- **Average Burst Pressure:** 30atm (Ø2.0-5.0mm); 26atm (Ø6.0-8.0mm)
- **Sheath Compatibility:** 4F (Ø2.0-4.0mm x 10-150mm; Ø4.5-6.0 x 10-60mm)
 5F (Ø4.5-6.0mm x 80-150mm; Ø7.0 x 10-150mm; Ø8.0-10.0 x 10-40mm)
 6F (Ø8.0 x 60-150mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
Catheter Working Length: 40cm				804-035010-040	10mm	3.5mm	40cm	804-050015-040	15mm	5.0mm	40cm
804-020010-040	10mm	2.0mm	40cm	804-035015-040	15mm	3.5mm	40cm	804-050020-040	20mm	5.0mm	40cm
804-020015-040	15mm	2.0mm	40cm	804-035020-040	20mm	3.5mm	40cm	804-050030-040	30mm	5.0mm	40cm
804-020020-040	20mm	2.0mm	40cm	804-035030-040	30mm	3.5mm	40cm	804-050040-040	40mm	5.0mm	40cm
804-020030-040	30mm	2.0mm	40cm	804-035040-040	40mm	3.5mm	40cm	804-050060-040	60mm	5.0mm	40cm
804-020040-040	40mm	2.0mm	40cm	804-035060-040	60mm	3.5mm	40cm	804-050080-040	80mm	5.0mm	40cm
804-020060-040	60mm	2.0mm	40cm	804-035080-040	80mm	3.5mm	40cm	804-050100-040	100mm	5.0mm	40cm
804-020080-040	80mm	2.0mm	40cm	804-035100-040	100mm	3.5mm	40cm	804-050150-040	150mm	5.0mm	40cm
804-020100-040	100mm	2.0mm	40cm	804-035150-040	150mm	3.5mm	40cm	804-055010-040	10mm	5.5mm	40cm
804-020150-040	150mm	2.0mm	40cm	804-040010-040	10mm	4.0mm	40cm	804-055015-040	15mm	5.5mm	40cm
804-025010-040	10mm	2.5mm	40cm	804-040015-040	15mm	4.0mm	40cm	804-055020-040	20mm	5.5mm	40cm
804-025015-040	15mm	2.5mm	40cm	804-040020-040	20mm	4.0mm	40cm	804-055030-040	30mm	5.5mm	40cm
804-025020-040	20mm	2.5mm	40cm	804-040030-040	30mm	4.0mm	40cm	804-055040-040	40mm	5.5mm	40cm
804-025030-040	30mm	2.5mm	40cm	804-040040-040	40mm	4.0mm	40cm	804-055060-040	60mm	5.5mm	40cm
804-025040-040	40mm	2.5mm	40cm	804-040060-040	60mm	4.0mm	40cm	804-055080-040	80mm	5.5mm	40cm
804-025060-040	60mm	2.5mm	40cm	804-040080-040	80mm	4.0mm	40cm	804-055100-040	100mm	5.5mm	40cm
804-025080-040	80mm	2.5mm	40cm	804-040100-040	100mm	4.0mm	40cm	804-055150-040	150mm	5.5mm	40cm
804-025100-040	100mm	2.5mm	40cm	804-040150-040	150mm	4.0mm	40cm	804-060010-040	10mm	6.0mm	40cm
804-025150-040	150mm	2.5mm	40cm	804-045010-040	10mm	4.5mm	40cm	804-060015-040	15mm	6.0mm	40cm
804-030010-040	10mm	3.0mm	40cm	804-045015-040	15mm	4.5mm	40cm	804-060020-040	20mm	6.0mm	40cm
804-030015-040	15mm	3.0mm	40cm	804-045020-040	20mm	4.5mm	40cm	804-060030-040	30mm	6.0mm	40cm
804-030020-040	20mm	3.0mm	40cm	804-045030-040	30mm	4.5mm	40cm	804-060040-040	40mm	6.0mm	40cm
804-030030-040	30mm	3.0mm	40cm	804-045040-040	40mm	4.5mm	40cm	804-060060-040	60mm	6.0mm	40cm
804-030040-040	40mm	3.0mm	40cm	804-045060-040	60mm	4.5mm	40cm	804-060080-040	80mm	6.0mm	40cm
804-030060-040	60mm	3.0mm	40cm	804-045080-040	80mm	4.5mm	40cm	804-060100-040	100mm	6.0mm	40cm
804-030080-040	80mm	3.0mm	40cm	804-045100-040	100mm	4.5mm	40cm	804-060150-040	150mm	6.0mm	40cm
804-030100-040	100mm	3.0mm	40cm	804-045150-040	150mm	4.5mm	40cm	804-070010-040	10mm	7.0mm	40cm
804-030150-040	150mm	3.0mm	40cm	804-050010-040	10mm	5.0mm	40cm	804-070015-040	15mm	7.0mm	40cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
804-070020-040	20mm	7.0mm	40cm	804-035020-070	20mm	3.5mm	70cm	804-060030-070	30mm	6.0mm	70cm
804-070030-040	30mm	7.0mm	40cm	804-035030-070	30mm	3.5mm	70cm	804-060040-070	40mm	6.0mm	70cm
804-070040-040	40mm	7.0mm	40cm	804-035040-070	40mm	3.5mm	70cm	804-060060-070	60mm	6.0mm	70cm
804-070060-040	60mm	7.0mm	40cm	804-035060-070	60mm	3.5mm	70cm	804-060080-070	80mm	6.0mm	70cm
804-070080-040	80mm	7.0mm	40cm	804-035080-070	80mm	3.5mm	70cm	804-060100-070	100mm	6.0mm	70cm
804-070100-040	100mm	7.0mm	40cm	804-035100-070	100mm	3.5mm	70cm	804-060150-070	150mm	6.0mm	70cm
804-070150-040	150mm	7.0mm	40cm	804-035150-070	150mm	3.5mm	70cm	804-070010-070	10mm	7.0mm	70cm
804-080010-040	10mm	8.0mm	40cm	804-040010-070	10mm	4.0mm	70cm	804-070015-070	15mm	7.0mm	70cm
804-080015-040	15mm	8.0mm	40cm	804-040015-070	15mm	4.0mm	70cm	804-070020-070	20mm	7.0mm	70cm
804-080020-040	20mm	8.0mm	40cm	804-040020-070	20mm	4.0mm	70cm	804-070030-070	30mm	7.0mm	70cm
804-080030-040	30mm	8.0mm	40cm	804-040030-070	30mm	4.0mm	70cm	804-070040-070	40mm	7.0mm	70cm
804-080040-040	40mm	8.0mm	40cm	804-040040-070	40mm	4.0mm	70cm	804-070060-070	60mm	7.0mm	70cm
804-080060-040	60mm	8.0mm	40cm	804-040060-070	60mm	4.0mm	70cm	804-070080-070	80mm	7.0mm	70cm
804-080080-040	80mm	8.0mm	40cm	804-040080-070	80mm	4.0mm	70cm	804-070100-070	100mm	7.0mm	70cm
804-080100-040	100mm	8.0mm	40cm	804-040100-070	100mm	4.0mm	70cm	804-070150-070	150mm	7.0mm	70cm
804-080150-040	150mm	8.0mm	40cm	804-040150-070	150mm	4.0mm	70cm	804-080010-070	10mm	8.0mm	70cm
Catheter Working Length: 70cm				804-045010-070	10mm	4.5mm	70cm	804-080015-070	15mm	8.0mm	70cm
804-020015-070	15mm	2.0mm	70cm	804-045015-070	15mm	4.5mm	70cm	804-080020-070	20mm	8.0mm	70cm
804-020020-070	20mm	2.0mm	70cm	804-045020-070	20mm	4.5mm	70cm	804-080030-070	30mm	8.0mm	70cm
804-020030-070	30mm	2.0mm	70cm	804-045030-070	30mm	4.5mm	70cm	804-080040-070	40mm	8.0mm	70cm
804-020040-070	40mm	2.0mm	70cm	804-045040-070	40mm	4.5mm	70cm	804-080060-070	60mm	8.0mm	70cm
804-020060-070	60mm	2.0mm	70cm	804-045060-070	60mm	4.5mm	70cm	804-080080-070	80mm	8.0mm	70cm
804-020080-070	80mm	2.0mm	70cm	804-045080-070	80mm	4.5mm	70cm	804-080100-070	100mm	8.0mm	70cm
804-020100-070	100mm	2.0mm	70cm	804-045100-070	100mm	4.5mm	70cm	804-080150-070	150mm	8.0mm	70cm
804-020150-070	150mm	2.0mm	70cm	804-045150-070	150mm	4.5mm	70cm	Catheter Working Length: 90cm			
804-025010-070	10mm	2.5mm	70cm	804-050010-070	10mm	5.0mm	70cm	804-020010-090	10mm	2.0mm	90cm
804-025015-070	15mm	2.5mm	70cm	804-050015-070	15mm	5.0mm	70cm	804-020015-090	15mm	2.0mm	90cm
804-025020-070	20mm	2.5mm	70cm	804-050020-070	20mm	5.0mm	70cm	804-020020-090	20mm	2.0mm	90cm
804-025030-070	30mm	2.5mm	70cm	804-050030-070	30mm	5.0mm	70cm	804-020030-090	30mm	2.0mm	90cm
804-025040-070	40mm	2.5mm	70cm	804-050040-070	40mm	5.0mm	70cm	804-020040-090	40mm	2.0mm	90cm
804-025060-070	60mm	2.5mm	70cm	804-050060-070	60mm	5.0mm	70cm	804-020060-090	60mm	2.0mm	90cm
804-025080-070	80mm	2.5mm	70cm	804-050080-070	80mm	5.0mm	70cm	804-020080-090	80mm	2.0mm	90cm
804-025100-070	100mm	2.5mm	70cm	804-050100-070	100mm	5.0mm	70cm	804-020100-090	100mm	2.0mm	90cm
804-025150-070	150mm	2.5mm	70cm	804-050150-070	150mm	5.0mm	70cm	804-020150-090	150mm	2.0mm	90cm
804-030010-070	10mm	3.0mm	70cm	804-055010-070	10mm	5.5mm	70cm	804-025010-090	10mm	2.5mm	90cm
804-030015-070	15mm	3.0mm	70cm	804-055015-070	15mm	5.5mm	70cm	804-025015-090	15mm	2.5mm	90cm
804-030020-070	20mm	3.0mm	70cm	804-055020-070	20mm	5.5mm	70cm	804-025020-090	20mm	2.5mm	90cm
804-030030-070	30mm	3.0mm	70cm	804-055030-070	30mm	5.5mm	70cm	804-025030-090	30mm	2.5mm	90cm
804-030040-070	40mm	3.0mm	70cm	804-055040-070	40mm	5.5mm	70cm	804-025040-090	40mm	2.5mm	90cm
804-030060-070	60mm	3.0mm	70cm	804-055060-070	60mm	5.5mm	70cm	804-025060-090	60mm	2.5mm	90cm
804-030080-070	80mm	3.0mm	70cm	804-055080-070	80mm	5.5mm	70cm	804-025080-090	80mm	2.5mm	90cm
804-030100-070	100mm	3.0mm	70cm	804-055100-070	100mm	5.5mm	70cm	804-025100-090	100mm	2.5mm	90cm
804-030150-070	150mm	3.0mm	70cm	804-055150-070	150mm	5.5mm	70cm	804-025150-090	150mm	2.5mm	90cm
804-035010-070	10mm	3.5mm	70cm	804-060010-070	10mm	6.0mm	70cm	804-030010-090	10mm	3.0mm	90cm
804-035015-070	15mm	3.5mm	70cm	804-060015-070	15mm	6.0mm	70cm	804-030015-090	15mm	3.0mm	90cm
				804-060020-070	20mm	6.0mm	70cm				

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
804-050030-150	30mm	5.0mm	150cm	804-055150-150	150mm	5.5mm	150cm	804-070040-150	40mm	7.0mm	150cm
804-050040-150	40mm	5.0mm	150cm	804-060010-150	10mm	6.0mm	150cm	804-070060-150	60mm	7.0mm	150cm
804-050060-150	60mm	5.0mm	150cm	804-060015-150	15mm	6.0mm	150cm	804-070080-150	80mm	7.0mm	150cm
804-050080-150	80mm	5.0mm	150cm	804-060020-150	20mm	6.0mm	150cm	804-070100-150	100mm	7.0mm	150cm
804-050100-150	100mm	5.0mm	150cm	804-060030-150	30mm	6.0mm	150cm	804-070150-150	150mm	7.0mm	150cm
804-050150-150	150mm	5.0mm	150cm	804-060040-150	40mm	6.0mm	150cm	804-080010-150	10mm	8.0mm	150cm
804-055010-150	10mm	5.5mm	150cm	804-060060-150	60mm	6.0mm	150cm	804-080015-150	15mm	8.0mm	150cm
804-055015-150	15mm	5.5mm	150cm	804-060080-150	80mm	6.0mm	150cm	804-080020-150	20mm	8.0mm	150cm
804-055020-150	20mm	5.5mm	150cm	804-060100-150	100mm	6.0mm	150cm	804-080030-150	30mm	8.0mm	150cm
804-055030-150	30mm	5.5mm	150cm	804-060150-150	150mm	6.0mm	150cm	804-080040-150	40mm	8.0mm	150cm
804-055040-150	40mm	5.5mm	150cm	804-070010-150	10mm	7.0mm	150cm	804-080060-150	60mm	8.0mm	150cm
804-055060-150	60mm	5.5mm	150cm	804-070015-150	15mm	7.0mm	150cm	804-080080-150	80mm	8.0mm	150cm
804-055080-150	80mm	5.5mm	150cm	804-070020-150	20mm	7.0mm	150cm	804-080100-150	100mm	8.0mm	150cm
804-055100-150	100mm	5.5mm	150cm	804-070030-150	30mm	7.0mm	150cm	804-080150-150	150mm	8.0mm	150cm

Zeeland

PTA Balloon Catheter

Semi-Compliant 0.014 OTW

- **Nominal Pressure:** 6atm
- **Rated Burst Pressure (RBP):** 14atm
- **Average Burst Pressure:** 20atm
- **Sheath Compatibility:** 4F (Ø1.5-5.0mm);
5F (Ø5.5-6.0mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
Catheter Working Length: 70cm				818-025020-070	20mm	2.5mm	70cm	818-035080-070	80mm	3.5mm	70cm
818-015005-070	5mm	1.5mm	70cm	818-025030-070	30mm	2.5mm	70cm	818-035100-070	100mm	3.5mm	70cm
818-015010-070	10mm	1.5mm	70cm	818-025040-070	40mm	2.5mm	70cm	818-035120-070	120mm	3.5mm	70cm
818-015015-070	15mm	1.5mm	70cm	818-025060-070	60mm	2.5mm	70cm	818-035150-070	150mm	3.5mm	70cm
818-015020-070	20mm	1.5mm	70cm	818-025080-070	80mm	2.5mm	70cm	818-035200-070	200mm	3.5mm	70cm
818-015030-070	30mm	1.5mm	70cm	818-025100-070	100mm	2.5mm	70cm	818-040005-070	5mm	4.0mm	70cm
818-015040-070	40mm	1.5mm	70cm	818-025120-070	120mm	2.5mm	70cm	818-040010-070	10mm	4.0mm	70cm
818-015060-070	60mm	1.5mm	70cm	818-025150-070	150mm	2.5mm	70cm	818-040015-070	15mm	4.0mm	70cm
818-015080-070	80mm	1.5mm	70cm	818-025200-070	200mm	2.5mm	70cm	818-040020-070	20mm	4.0mm	70cm
818-015100-070	100mm	1.5mm	70cm	818-030005-070	5mm	3.0mm	70cm	818-040030-070	30mm	4.0mm	70cm
818-015120-070	120mm	1.5mm	70cm	818-030010-070	10mm	3.0mm	70cm	818-040040-070	40mm	4.0mm	70cm
818-015150-070	150mm	1.5mm	70cm	818-030015-070	15mm	3.0mm	70cm	818-040060-070	60mm	4.0mm	70cm
818-015200-070	200mm	1.5mm	70cm	818-030020-070	20mm	3.0mm	70cm	818-040080-070	80mm	4.0mm	70cm
818-020005-070	5mm	2.0mm	70cm	818-030030-070	30mm	3.0mm	70cm	818-040100-070	100mm	4.0mm	70cm
818-020010-070	10mm	2.0mm	70cm	818-030040-070	40mm	3.0mm	70cm	818-040120-070	120mm	4.0mm	70cm
818-020015-070	15mm	2.0mm	70cm	818-030060-070	60mm	3.0mm	70cm	818-040150-070	150mm	4.0mm	70cm
818-020020-070	20mm	2.0mm	70cm	818-030080-070	80mm	3.0mm	70cm	818-040200-070	200mm	4.0mm	70cm
818-020030-070	30mm	2.0mm	70cm	818-030100-070	100mm	3.0mm	70cm	818-045005-070	5mm	4.5mm	70cm
818-020040-070	40mm	2.0mm	70cm	818-030120-070	120mm	3.0mm	70cm	818-045010-070	10mm	4.5mm	70cm
818-020060-070	60mm	2.0mm	70cm	818-030150-070	150mm	3.0mm	70cm	818-045015-070	15mm	4.5mm	70cm
818-020080-070	80mm	2.0mm	70cm	818-030200-070	200mm	3.0mm	70cm	818-045020-070	20mm	4.5mm	70cm
818-020100-070	100mm	2.0mm	70cm	818-035005-070	5mm	3.5mm	70cm	818-045030-070	30mm	4.5mm	70cm
818-020120-070	120mm	2.0mm	70cm	818-035010-070	10mm	3.5mm	70cm	818-045040-070	40mm	4.5mm	70cm
818-020150-070	150mm	2.0mm	70cm	818-035015-070	15mm	3.5mm	70cm	818-045060-070	60mm	4.5mm	70cm
818-020200-070	200mm	2.0mm	70cm	818-035020-070	20mm	3.5mm	70cm	818-045080-070	80mm	4.5mm	70cm
818-025005-070	5mm	2.5mm	70cm	818-035030-070	30mm	3.5mm	70cm	818-045100-070	100mm	4.5mm	70cm
818-025010-070	10mm	2.5mm	70cm	818-035040-070	40mm	3.5mm	70cm	818-045120-070	120mm	4.5mm	70cm
818-025015-070	15mm	2.5mm	70cm	818-035060-070	60mm	3.5mm	70cm	818-045150-070	150mm	4.5mm	70cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
818-055020-090	20mm	5.5mm	90cm	818-025010-150	10mm	2.5mm	150cm	818-045010-150	10mm	4.5mm	150cm
818-055030-090	30mm	5.5mm	90cm	818-025015-150	15mm	2.5mm	150cm	818-045015-150	15mm	4.5mm	150cm
818-055040-090	40mm	5.5mm	90cm	818-025020-150	20mm	2.5mm	150cm	818-045020-150	20mm	4.5mm	150cm
818-055060-090	60mm	5.5mm	90cm	818-025030-150	30mm	2.5mm	150cm	818-045030-150	30mm	4.5mm	150cm
818-055080-090	80mm	5.5mm	90cm	818-025040-150	40mm	2.5mm	150cm	818-045040-150	40mm	4.5mm	150cm
818-055100-090	100mm	5.5mm	90cm	818-025060-150	60mm	2.5mm	150cm	818-045060-150	60mm	4.5mm	150cm
818-055120-090	120mm	5.5mm	90cm	818-025080-150	80mm	2.5mm	150cm	818-045080-150	80mm	4.5mm	150cm
818-055150-090	150mm	5.5mm	90cm	818-025100-150	100mm	2.5mm	150cm	818-045100-150	100mm	4.5mm	150cm
818-055200-090	200mm	5.5mm	90cm	818-025120-150	120mm	2.5mm	150cm	818-045120-150	120mm	4.5mm	150cm
818-060005-090	5mm	6.0mm	90cm	818-025150-150	150mm	2.5mm	150cm	818-045150-150	150mm	4.5mm	150cm
818-060010-090	10mm	6.0mm	90cm	818-025200-150	200mm	2.5mm	150cm	818-045200-150	200mm	4.5mm	150cm
818-060015-090	15mm	6.0mm	90cm	818-030005-150	5mm	3.0mm	150cm	818-050005-150	5mm	5.0mm	150cm
818-060020-090	20mm	6.0mm	90cm	818-030010-150	10mm	3.0mm	150cm	818-050010-150	10mm	5.0mm	150cm
818-060030-090	30mm	6.0mm	90cm	818-030015-150	15mm	3.0mm	150cm	818-050015-150	15mm	5.0mm	150cm
818-060040-090	40mm	6.0mm	90cm	818-030020-150	20mm	3.0mm	150cm	818-050020-150	20mm	5.0mm	150cm
818-060060-090	60mm	6.0mm	90cm	818-030030-150	30mm	3.0mm	150cm	818-050030-150	30mm	5.0mm	150cm
818-060080-090	80mm	6.0mm	90cm	818-030040-150	40mm	3.0mm	150cm	818-050040-150	40mm	5.0mm	150cm
818-060100-090	100mm	6.0mm	90cm	818-030060-150	60mm	3.0mm	150cm	818-050060-150	60mm	5.0mm	150cm
818-060120-090	120mm	6.0mm	90cm	818-030080-150	80mm	3.0mm	150cm	818-050080-150	80mm	5.0mm	150cm
818-060150-090	150mm	6.0mm	90cm	818-030100-150	100mm	3.0mm	150cm	818-050100-150	100mm	5.0mm	150cm
818-060200-090	200mm	6.0mm	90cm	818-030120-150	120mm	3.0mm	150cm	818-050120-150	120mm	5.0mm	150cm
				818-030150-150	150mm	3.0mm	150cm	818-050150-150	150mm	5.0mm	150cm
				818-030200-150	200mm	3.0mm	150cm	818-050200-150	200mm	5.0mm	150cm
				818-035005-150	5mm	3.5mm	150cm	818-055005-150	5mm	5.5mm	150cm
				818-035010-150	10mm	3.5mm	150cm	818-055010-150	10mm	5.5mm	150cm
				818-035015-150	15mm	3.5mm	150cm	818-055015-150	15mm	5.5mm	150cm
				818-035020-150	20mm	3.5mm	150cm	818-055020-150	20mm	5.5mm	150cm
				818-035030-150	30mm	3.5mm	150cm	818-055030-150	30mm	5.5mm	150cm
				818-035040-150	40mm	3.5mm	150cm	818-055040-150	40mm	5.5mm	150cm
				818-035060-150	60mm	3.5mm	150cm	818-055060-150	60mm	5.5mm	150cm
				818-035080-150	80mm	3.5mm	150cm	818-055080-150	80mm	5.5mm	150cm
				818-035100-150	100mm	3.5mm	150cm	818-055100-150	100mm	5.5mm	150cm
				818-035120-150	120mm	3.5mm	150cm	818-055120-150	120mm	5.5mm	150cm
				818-035150-150	150mm	3.5mm	150cm	818-055150-150	150mm	5.5mm	150cm
				818-035200-150	200mm	3.5mm	150cm	818-055200-150	200mm	5.5mm	150cm
				818-040005-150	5mm	4.0mm	150cm	818-060005-150	5mm	6.0mm	150cm
				818-040010-150	10mm	4.0mm	150cm	818-060010-150	10mm	6.0mm	150cm
				818-040015-150	15mm	4.0mm	150cm	818-060015-150	15mm	6.0mm	150cm
				818-040020-150	20mm	4.0mm	150cm	818-060020-150	20mm	6.0mm	150cm
				818-040030-150	30mm	4.0mm	150cm	818-060030-150	30mm	6.0mm	150cm
				818-040040-150	40mm	4.0mm	150cm	818-060040-150	40mm	6.0mm	150cm
				818-040060-150	60mm	4.0mm	150cm	818-060060-150	60mm	6.0mm	150cm
				818-040080-150	80mm	4.0mm	150cm	818-060080-150	80mm	6.0mm	150cm
				818-040100-150	100mm	4.0mm	150cm	818-060100-150	100mm	6.0mm	150cm
				818-040120-150	120mm	4.0mm	150cm	818-060120-150	120mm	6.0mm	150cm
				818-040150-150	150mm	4.0mm	150cm	818-060150-150	150mm	6.0mm	150cm
				818-040200-150	200mm	4.0mm	150cm	818-060200-150	200mm	6.0mm	150cm
				818-045005-150	5mm	4.5mm	150cm				

Catheter Working Length: 150cm





Rialto NC

PTA Balloon Catheter

Non-Compliant 0.014 OTW

- **Nominal Pressure:** 12atm
- **Rated Burst Pressure (RBP):** 22atm (Ø2.0-4.0mm); 20atm (Ø5.0-6.0mm); 18atm (Ø7.0-8.0mm)
- **Average Burst Pressure:** 30atm (Ø2.0-5.0mm); 26atm (Ø6.0-8.0mm)
- **Sheath Compatibility:** 4F (Ø2.0-5.0mm); 5F (Ø5.5-6.0mm); 6F (Ø7.0-8.0mm)

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
816-020010-070	10mm	2.0mm	70cm	816-035015-070	15mm	3.5mm	70cm	816-050020-070	20mm	5.0mm	70cm
816-020015-070	15mm	2.0mm	70cm	816-035020-070	20mm	3.5mm	70cm	816-050030-070	30mm	5.0mm	70cm
816-020020-070	20mm	2.0mm	70cm	816-035030-070	30mm	3.5mm	70cm	816-050040-070	40mm	5.0mm	70cm
816-020030-070	30mm	2.0mm	70cm	816-035040-070	40mm	3.5mm	70cm	816-050060-070	60mm	5.0mm	70cm
816-020040-070	40mm	2.0mm	70cm	816-035060-070	60mm	3.5mm	70cm	816-050080-070	80mm	5.0mm	70cm
816-020060-070	60mm	2.0mm	70cm	816-035080-070	80mm	3.5mm	70cm	816-050100-070	100mm	5.0mm	70cm
816-020080-070	80mm	2.0mm	70cm	816-035100-070	100mm	3.5mm	70cm	816-050150-070	150mm	5.0mm	70cm
816-020100-070	100mm	2.0mm	70cm	816-035150-070	150mm	3.5mm	70cm	816-055010-070	10mm	5.5mm	70cm
816-020150-070	150mm	2.0mm	70cm	816-040010-070	10mm	4.0mm	70cm	816-055015-070	15mm	5.5mm	70cm
816-025010-070	10mm	2.5mm	70cm	816-040015-070	15mm	4.0mm	70cm	816-055020-070	20mm	5.5mm	70cm
816-025015-070	15mm	2.5mm	70cm	816-040020-070	20mm	4.0mm	70cm	816-055030-070	30mm	5.5mm	70cm
816-025020-070	20mm	2.5mm	70cm	816-040030-070	30mm	4.0mm	70cm	816-055040-070	40mm	5.5mm	70cm
816-025030-070	30mm	2.5mm	70cm	816-040040-070	40mm	4.0mm	70cm	816-055060-070	60mm	5.5mm	70cm
816-025040-070	40mm	2.5mm	70cm	816-040060-070	60mm	4.0mm	70cm	816-055080-070	80mm	5.5mm	70cm
816-025060-070	60mm	2.5mm	70cm	816-040080-070	80mm	4.0mm	70cm	816-055100-070	100mm	5.5mm	70cm
816-025080-070	80mm	2.5mm	70cm	816-040100-070	100mm	4.0mm	70cm	816-055150-070	150mm	5.5mm	70cm
816-025100-070	100mm	2.5mm	70cm	816-040150-070	150mm	4.0mm	70cm	816-060010-070	10mm	6.0mm	70cm
816-025150-070	150mm	2.5mm	70cm	816-045010-070	10mm	4.5mm	70cm	816-060015-070	15mm	6.0mm	70cm
816-030010-070	10mm	3.0mm	70cm	816-045015-070	15mm	4.5mm	70cm	816-060020-070	20mm	6.0mm	70cm
816-030015-070	15mm	3.0mm	70cm	816-045020-070	20mm	4.5mm	70cm	816-060030-070	30mm	6.0mm	70cm
816-030020-070	20mm	3.0mm	70cm	816-045030-070	30mm	4.5mm	70cm	816-060040-070	40mm	6.0mm	70cm
816-030030-070	30mm	3.0mm	70cm	816-045040-070	40mm	4.5mm	70cm	816-060060-070	60mm	6.0mm	70cm
816-030040-070	40mm	3.0mm	70cm	816-045060-070	60mm	4.5mm	70cm	816-060080-070	80mm	6.0mm	70cm
816-030060-070	60mm	3.0mm	70cm	816-045080-070	80mm	4.5mm	70cm	816-060100-070	100mm	6.0mm	70cm
816-030080-070	80mm	3.0mm	70cm	816-045100-070	100mm	4.5mm	70cm	816-060150-070	150mm	6.0mm	70cm
816-030100-070	100mm	3.0mm	70cm	816-045150-070	150mm	4.5mm	70cm	816-070010-070	10mm	7.0mm	70cm
816-030150-070	150mm	3.0mm	70cm	816-050010-070	10mm	5.0mm	70cm	816-070015-070	15mm	7.0mm	70cm
816-035010-070	10mm	3.5mm	70cm	816-050015-070	15mm	5.0mm	70cm	816-070020-070	20mm	7.0mm	70cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
816-070030-070	30mm	7.0mm	70cm	816-040020-090	20mm	4.0mm	90cm	816-070030-090	30mm	7.0mm	90cm
816-070040-070	40mm	7.0mm	70cm	816-040030-090	30mm	4.0mm	90cm	816-070040-090	40mm	7.0mm	90cm
816-070060-070	60mm	7.0mm	70cm	816-040040-090	40mm	4.0mm	90cm	816-070060-090	60mm	7.0mm	90cm
816-070080-070	80mm	7.0mm	70cm	816-040060-090	60mm	4.0mm	90cm	816-070080-090	80mm	7.0mm	90cm
816-070100-070	100mm	7.0mm	70cm	816-040080-090	80mm	4.0mm	90cm	816-070100-090	100mm	7.0mm	90cm
816-070150-070	150mm	7.0mm	70cm	816-040100-090	100mm	4.0mm	90cm	816-070150-090	150mm	7.0mm	90cm
				816-040150-090	150mm	4.0mm	90cm				
Catheter Working Length: 90cm				816-045010-090	10mm	4.5mm	90cm	Catheter Working Length: 150cm			
816-020010-090	10mm	2.0mm	90cm	816-045015-090	15mm	4.5mm	90cm	816-020010-150	10mm	2.0mm	150cm
816-020015-090	15mm	2.0mm	90cm	816-045020-090	20mm	4.5mm	90cm	816-020015-150	15mm	2.0mm	150cm
816-020020-090	20mm	2.0mm	90cm	816-045030-090	30mm	4.5mm	90cm	816-020020-150	20mm	2.0mm	150cm
816-020030-090	30mm	2.0mm	90cm	816-045040-090	40mm	4.5mm	90cm	816-020030-150	30mm	2.0mm	150cm
816-020040-090	40mm	2.0mm	90cm	816-045060-090	60mm	4.5mm	90cm	816-020040-150	40mm	2.0mm	150cm
816-020060-090	60mm	2.0mm	90cm	816-045080-090	80mm	4.5mm	90cm	816-020060-150	60mm	2.0mm	150cm
816-020080-090	80mm	2.0mm	90cm	816-045100-090	100mm	4.5mm	90cm	816-020080-150	80mm	2.0mm	150cm
816-020100-090	100mm	2.0mm	90cm	816-045150-090	150mm	4.5mm	90cm	816-020100-150	100mm	2.0mm	150cm
816-020150-090	150mm	2.0mm	90cm	816-050010-090	10mm	5.0mm	90cm	816-020150-150	150mm	2.0mm	150cm
816-025010-090	10mm	2.5mm	90cm	816-050015-090	15mm	5.0mm	90cm	816-025010-150	10mm	2.5mm	150cm
816-025015-090	15mm	2.5mm	90cm	816-050020-090	20mm	5.0mm	90cm	816-025015-150	15mm	2.5mm	150cm
816-025020-090	20mm	2.5mm	90cm	816-050030-090	30mm	5.0mm	90cm	816-025020-150	20mm	2.5mm	150cm
816-025030-090	30mm	2.5mm	90cm	816-050040-090	40mm	5.0mm	90cm	816-025030-150	30mm	2.5mm	150cm
816-025040-090	40mm	2.5mm	90cm	816-050060-090	60mm	5.0mm	90cm	816-025040-150	40mm	2.5mm	150cm
816-025060-090	60mm	2.5mm	90cm	816-050080-090	80mm	5.0mm	90cm	816-025060-150	60mm	2.5mm	150cm
816-025080-090	80mm	2.5mm	90cm	816-050100-090	100mm	5.0mm	90cm	816-025080-150	80mm	2.5mm	150cm
816-025100-090	100mm	2.5mm	90cm	816-050150-090	150mm	5.0mm	90cm	816-025100-150	100mm	2.5mm	150cm
816-025150-090	150mm	2.5mm	90cm	816-055010-090	10mm	5.5mm	90cm	816-025150-150	150mm	2.5mm	150cm
816-030010-090	10mm	3.0mm	90cm	816-055015-090	15mm	5.5mm	90cm	816-030010-150	10mm	3.0mm	150cm
816-030015-090	15mm	3.0mm	90cm	816-055020-090	20mm	5.5mm	90cm	816-030015-150	15mm	3.0mm	150cm
816-030020-090	20mm	3.0mm	90cm	816-055030-090	30mm	5.5mm	90cm	816-030020-150	20mm	3.0mm	150cm
816-030030-090	30mm	3.0mm	90cm	816-055040-090	40mm	5.5mm	90cm	816-030030-150	30mm	3.0mm	150cm
816-030040-090	40mm	3.0mm	90cm	816-055060-090	60mm	5.5mm	90cm	816-030040-150	40mm	3.0mm	150cm
816-030060-090	60mm	3.0mm	90cm	816-055080-090	80mm	5.5mm	90cm	816-030060-150	60mm	3.0mm	150cm
816-030080-090	80mm	3.0mm	90cm	816-055100-090	100mm	5.5mm	90cm	816-030080-150	80mm	3.0mm	150cm
816-030100-090	100mm	3.0mm	90cm	816-055150-090	150mm	5.5mm	90cm	816-030100-150	100mm	3.0mm	150cm
816-030150-090	150mm	3.0mm	90cm	816-060010-090	10mm	6.0mm	90cm	816-030150-150	150mm	3.0mm	150cm
816-035010-090	10mm	3.5mm	90cm	816-060015-090	15mm	6.0mm	90cm	816-035010-150	10mm	3.5mm	150cm
816-035015-090	15mm	3.5mm	90cm	816-060020-090	20mm	6.0mm	90cm	816-035015-150	15mm	3.5mm	150cm
816-035020-090	20mm	3.5mm	90cm	816-060030-090	30mm	6.0mm	90cm	816-035020-150	20mm	3.5mm	150cm
816-035030-090	30mm	3.5mm	90cm	816-060040-090	40mm	6.0mm	90cm	816-035030-150	30mm	3.5mm	150cm
816-035040-090	40mm	3.5mm	90cm	816-060060-090	60mm	6.0mm	90cm	816-035040-150	40mm	3.5mm	150cm
816-035060-090	60mm	3.5mm	90cm	816-060080-090	80mm	6.0mm	90cm	816-035060-150	60mm	3.5mm	150cm
816-035080-090	80mm	3.5mm	90cm	816-060100-090	100mm	6.0mm	90cm	816-035080-150	80mm	3.5mm	150cm
816-035100-090	100mm	3.5mm	90cm	816-060150-090	150mm	6.0mm	90cm	816-035100-150	100mm	3.5mm	150cm
816-035150-090	150mm	3.5mm	90cm	816-070010-090	10mm	7.0mm	90cm	816-035150-150	150mm	3.5mm	150cm
816-040010-090	10mm	4.0mm	90cm	816-070015-090	15mm	7.0mm	90cm	816-040010-150	10mm	4.0mm	150cm
816-040015-090	15mm	4.0mm	90cm	816-070020-090	20mm	7.0mm	90cm	816-040015-150	15mm	4.0mm	150cm

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
816-040020-150	20mm	4.0mm	150cm	816-050020-150	20mm	5.0mm	150cm	816-060020-150	20mm	6.0mm	150cm
816-040030-150	30mm	4.0mm	150cm	816-050030-150	30mm	5.0mm	150cm	816-060030-150	30mm	6.0mm	150cm
816-040040-150	40mm	4.0mm	150cm	816-050040-150	40mm	5.0mm	150cm	816-060040-150	40mm	6.0mm	150cm
816-040060-150	60mm	4.0mm	150cm	816-050060-150	60mm	5.0mm	150cm	816-060060-150	60mm	6.0mm	150cm
816-040080-150	80mm	4.0mm	150cm	816-050080-150	80mm	5.0mm	150cm	816-060080-150	80mm	6.0mm	150cm
816-040100-150	100mm	4.0mm	150cm	816-050100-150	100mm	5.0mm	150cm	816-060100-150	100mm	6.0mm	150cm
816-040150-150	150mm	4.0mm	150cm	816-050150-150	150mm	5.0mm	150cm	816-060150-150	150mm	6.0mm	150cm
816-045010-150	10mm	4.5mm	150cm	816-055010-150	10mm	5.5mm	150cm	816-070010-150	10mm	7.0mm	150cm
816-045015-150	15mm	4.5mm	150cm	816-055015-150	15mm	5.5mm	150cm	816-070015-150	15mm	7.0mm	150cm
816-045020-150	20mm	4.5mm	150cm	816-055020-150	20mm	5.5mm	150cm	816-070020-150	20mm	7.0mm	150cm
816-045030-150	30mm	4.5mm	150cm	816-055030-150	30mm	5.5mm	150cm	816-070030-150	30mm	7.0mm	150cm
816-045040-150	40mm	4.5mm	150cm	816-055040-150	40mm	5.5mm	150cm	816-070040-150	40mm	7.0mm	150cm
816-045060-150	60mm	4.5mm	150cm	816-055060-150	60mm	5.5mm	150cm	816-070060-150	60mm	7.0mm	150cm
816-045080-150	80mm	4.5mm	150cm	816-055080-150	80mm	5.5mm	150cm	816-070080-150	80mm	7.0mm	150cm
816-045100-150	100mm	4.5mm	150cm	816-055100-150	100mm	5.5mm	150cm	816-070100-150	100mm	7.0mm	150cm
816-045150-150	150mm	4.5mm	150cm	816-055150-150	150mm	5.5mm	150cm	816-070150-150	150mm	7.0mm	150cm
816-050010-150	10mm	5.0mm	150cm	816-060010-150	10mm	6.0mm	150cm				
816-050015-150	15mm	5.0mm	150cm	816-060015-150	15mm	6.0mm	150cm				





Helix

PTCA Balloon Catheter Semi-Compliant

- **Nominal Pressure:** 6atm
- **Rated Burst Pressure (RBP):** 14atm
- **Average Burst Pressure:** 20atm
- **Sheath Compatibility:** 5F compatible / 6F compatible for kissing

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
801-1005	5mm	1.0mm	140cm	801-1715	15mm	1.75mm	140cm	801-2705	5mm	2.75mm	140cm
801-1008	8mm	1.0mm	140cm*	801-1720	20mm	1.75mm	140cm	801-2708	8mm	2.75mm	140cm*
801-1010	10mm	1.0mm	140cm	801-1725	25mm	1.75mm	140cm*	801-2710	10mm	2.75mm	140cm
801-1012	12mm	1.0mm	140cm*	801-1730	30mm	1.75mm	140cm	801-2712	12mm	2.75mm	140cm*
801-1015	15mm	1.0mm	140cm	801-2005	5mm	2.0mm	140cm	801-2715	15mm	2.75mm	140cm
801-1020	20mm	1.0mm	140cm	801-2008	8mm	2.0mm	140cm*	801-2720	20mm	2.75mm	140cm
801-1025	25mm	1.0mm	140cm*	801-2010	10mm	2.0mm	140cm	801-2725	25mm	2.75mm	140cm*
801-1030	30mm	1.0mm	140cm	801-2012	12mm	2.0mm	140cm*	801-2730	30mm	2.75mm	140cm
801-1205	5mm	1.25mm	140cm	801-2015	15mm	2.0mm	140cm	801-3005	5mm	3.0mm	140cm
801-1208	8mm	1.25mm	140cm*	801-2020	20mm	2.0mm	140cm	801-3008	8mm	3.0mm	140cm*
801-1210	10mm	1.25mm	140cm	801-2025	25mm	2.0mm	140cm*	801-3010	10mm	3.0mm	140cm
801-1212	12mm	1.25mm	140cm*	801-2030	30mm	2.0mm	140cm	801-3012	12mm	3.0mm	140cm*
801-1215	15mm	1.25mm	140cm	801-2205	5mm	2.25mm	140cm	801-3015	15mm	3.0mm	140cm
801-1220	20mm	1.25mm	140cm	801-2208	8mm	2.25mm	140cm*	801-3020	20mm	3.0mm	140cm
801-1225	25mm	1.25mm	140cm*	801-2210	10mm	2.25mm	140cm	801-3025	25mm	3.0mm	140cm*
801-1230	30mm	1.25mm	140cm	801-2212	12mm	2.25mm	140cm*	801-3030	30mm	3.0mm	140cm
801-1505	5mm	1.5mm	140cm	801-2215	15mm	2.25mm	140cm	801-3205	5mm	3.25mm	140cm
801-1508	8mm	1.5mm	140cm*	801-2220	20mm	2.25mm	140cm	801-3208	8mm	3.25mm	140cm*
801-1510	10mm	1.5mm	140cm	801-2225	25mm	2.25mm	140cm*	801-3210	10mm	3.25mm	140cm
801-1512	12mm	1.5mm	140cm*	801-2230	30mm	2.25mm	140cm	801-3212	12mm	3.25mm	140cm*
801-1515	15mm	1.5mm	140cm	801-2505	5mm	2.5mm	140cm	801-3215	15mm	3.25mm	140cm
801-1520	20mm	1.5mm	140cm	801-2508	8mm	2.5mm	140cm*	801-3220	20mm	3.25mm	140cm
801-1525	25mm	1.5mm	140cm*	801-2510	10mm	2.5mm	140cm	801-3225	25mm	3.25mm	140cm*
801-1530	30mm	1.5mm	140cm	801-2512	12mm	2.5mm	140cm*	801-3230	30mm	3.25mm	140cm
801-1705	5mm	1.75mm	140cm	801-2515	15mm	2.5mm	140cm	801-3505	5mm	3.5mm	140cm
801-1708	8mm	1.75mm	140cm*	801-2520	20mm	2.5mm	140cm	801-3508	8mm	3.5mm	140cm*
801-1710	10mm	1.75mm	140cm	801-2525	25mm	2.5mm	140cm*	801-3510	10mm	3.5mm	140cm
801-1712	12mm	1.75mm	140cm*	801-2530	30mm	2.5mm	140cm	801-3512	12mm	3.5mm	140cm*

Note: * available after CE approval

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
801-3515	15mm	3.5mm	140cm	801-3715	15mm	3.75mm	140cm	801-4015	15mm	4.0mm	140cm
801-3520	20mm	3.5mm	140cm	801-3720	20mm	3.75mm	140cm	801-4020	20mm	4.0mm	140cm
801-3525	25mm	3.5mm	140cm*	801-3725	25mm	3.75mm	140cm*	801-4025	25mm	4.0mm	140cm*
801-3530	30mm	3.5mm	140cm	801-3730	30mm	3.75mm	140cm	801-4030	30mm	4.0mm	140cm
801-3705	5mm	3.75mm	140cm	801-4005	5mm	4.0mm	140cm				
801-3708	8mm	3.75mm	140cm*	801-4008	8mm	4.0mm	140cm*				
801-3710	10mm	3.75mm	140cm	801-4010	10mm	4.0mm	140cm				
801-3712	12mm	3.75mm	140cm*	801-4012	12mm	4.0mm	140cm*				

Note: * available after CE approval



Vecchio NC

PTCA Balloon Catheter Non-Compliant

- **Nominal Pressure:** 12atm
- **Rated Burst Pressure (RBP):** 22atm (Ø2.0-4.0mm); 20atm (Ø4.5-5.0mm)
- **Average Burst Pressure:** 30atm
- **Sheath Compatibility:** 5F compatible / 6F compatible for kissing

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
802-2006	6mm	2.0mm	140cm*	802-2525	25mm	2.5mm	140cm*	802-3220	20mm	3.25mm	140cm*
802-2008	8mm	2.0mm	140cm	802-2530	30mm	2.5mm	140cm*	802-3222	22mm	3.25mm	140cm*
802-2010	10mm	2.0mm	140cm	802-2706	6mm	2.75mm	140cm*	802-3225	25mm	3.25mm	140cm*
802-2012	12mm	2.0mm	140cm	802-2708	8mm	2.75mm	140cm	802-3230	30mm	3.25mm	140cm*
802-2015	15mm	2.0mm	140cm	802-2710	10mm	2.75mm	140cm	802-3506	6mm	3.5mm	140cm*
802-2018	18mm	2.0mm	140cm	802-2712	12mm	2.75mm	140cm	802-3508	8mm	3.5mm	140cm
802-2020	20mm	2.0mm	140cm*	802-2715	15mm	2.75mm	140cm	802-3510	10mm	3.5mm	140cm
802-2022	22mm	2.0mm	140cm*	802-2718	18mm	2.75mm	140cm	802-3512	12mm	3.5mm	140cm
802-2025	25mm	2.0mm	140cm*	802-2720	20mm	2.75mm	140cm*	802-3515	15mm	3.5mm	140cm
802-2030	30mm	2.0mm	140cm*	802-2722	22mm	2.75mm	140cm*	802-3518	18mm	3.5mm	140cm
802-2206	6mm	2.25mm	140cm*	802-2725	25mm	2.75mm	140cm*	802-3520	20mm	3.5mm	140cm*
802-2208	8mm	2.25mm	140cm	802-2730	30mm	2.75mm	140cm*	802-3522	22mm	3.5mm	140cm*
802-2210	10mm	2.25mm	140cm	802-3006	6mm	3.0mm	140cm*	802-3525	25mm	3.5mm	140cm*
802-2212	12mm	2.25mm	140cm	802-3008	8mm	3.0mm	140cm	802-3530	30mm	3.5mm	140cm*
802-2215	15mm	2.25mm	140cm	802-3010	10mm	3.0mm	140cm	802-3706	6mm	3.75mm	140cm*
802-2218	18mm	2.25mm	140cm	802-3012	12mm	3.0mm	140cm	802-3708	8mm	3.75mm	140cm
802-2220	20mm	2.25mm	140cm*	802-3015	15mm	3.0mm	140cm	802-3710	10mm	3.75mm	140cm
802-2222	22mm	2.25mm	140cm*	802-3018	18mm	3.0mm	140cm	802-3712	12mm	3.75mm	140cm
802-2225	25mm	2.25mm	140cm*	802-3020	20mm	3.0mm	140cm*	802-3715	15mm	3.75mm	140cm
802-2230	30mm	2.25mm	140cm*	802-3022	22mm	3.0mm	140cm*	802-3718	18mm	3.75mm	140cm
802-2506	6mm	2.5mm	140cm*	802-3025	25mm	3.0mm	140cm*	802-3720	20mm	3.75mm	140cm*
802-2508	8mm	2.5mm	140cm	802-3030	30mm	3.0mm	140cm*	802-3722	22mm	3.75mm	140cm*
802-2510	10mm	2.5mm	140cm	802-3206	6mm	3.25mm	140cm*	802-3725	25mm	3.75mm	140cm*
802-2512	12mm	2.5mm	140cm	802-3208	8mm	3.25mm	140cm	802-3730	30mm	3.75mm	140cm*
802-2515	15mm	2.5mm	140cm	802-3210	10mm	3.25mm	140cm	802-4006	6mm	4.0mm	140cm*
802-2518	18mm	2.5mm	140cm	802-3212	12mm	3.25mm	140cm	802-4008	8mm	4.0mm	140cm
802-2520	20mm	2.5mm	140cm*	802-3215	15mm	3.25mm	140cm	802-4010	10mm	4.0mm	140cm
802-2522	22mm	2.5mm	140cm*	802-3218	18mm	3.25mm	140cm	802-4012	12mm	4.0mm	140cm

Note: * available after CE approval

Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length	Catalog No.	balloon Length	Balloon OD	Catheter Length
802-4015	15mm	4.0mm	140cm	802-4515	15mm	4.5mm	140cm	802-5015	15mm	5.0mm	140cm
802-4018	18mm	4.0mm	140cm	802-4518	18mm	4.5mm	140cm	802-5018	18mm	5.0mm	140cm
802-4020	20mm	4.0mm	140cm*	802-4520	20mm	4.5mm	140cm*	802-5020	20mm	5.0mm	140cm*
802-4022	22mm	4.0mm	140cm*	802-4522	22mm	4.5mm	140cm*	802-5022	22mm	5.0mm	140cm*
802-4025	25mm	4.0mm	140cm*	802-4525	25mm	4.5mm	140cm*	802-5025	25mm	5.0mm	140cm*
802-4030	30mm	4.0mm	140cm*	802-4530	30mm	4.5mm	140cm*	802-5030	30mm	5.0mm	140cm*
802-4506	6mm	4.5mm	140cm*	802-5006	6mm	5.0mm	140cm*				
802-4508	8mm	4.5mm	140cm	802-5008	8mm	5.0mm	140cm				
802-4510	10mm	4.5mm	140cm	802-5010	10mm	5.0mm	140cm				
802-4512	12mm	4.5mm	140cm	802-5012	12mm	5.0mm	140cm				

Note: * available after CE approval



Address:

cNovate Medical BV, Hogebrinkerweg 31, 3871 KM Hoevelaken, The Netherlands
t. +31 (0)850 140404, e. cs@cNovate.eu